


# ophtha futur<sup>®</sup> gas



High purity

Adjustable gas concentration

Environmentally conscious system

24 months shelf life

# ophthafutur<sup>®</sup> gas

## technical data-sheet

Product Name	ophthafutur <sup>®</sup> sf6	ophthafutur <sup>®</sup> c2f6	ophthafutur <sup>®</sup> c3f8
Active compound	sulfur hexafluoride	hexafluoroethane	octafluoropropane
<b>Composition</b>			
Sum Formula	SF <sub>6</sub>	C <sub>2</sub> F <sub>6</sub>	C <sub>3</sub> F <sub>8</sub>
Purity	> 4.5	> 5.0	> 4.0
Mass fraction	> 99.995%	> 99.999%	> 99.99%
<b>Non expansive gas concentration</b>	<b>20%</b>	<b>16%</b>	<b>12%</b>
<b>Potential tamponade time</b>	6 days	15 days	30 days
<b>Retention time</b>	10-14 days	30-35 days	55-65 days
<b>Properties</b>			
Pure Gas Volume	15 ml	12 ml	9 ml
Total Syringe Volume - Mixing Device	60 ml	60 ml	60 ml
<b>Shelf life</b>	24 months	24 months	24 months
<b>Presentation</b>			
Container	Glass Syringe	Glass Syringe	Glass Syringe
Mixing Device	Plastic Syringe with Filter	Plastic Syringe with Filter	Plastic Syringe with Filter
<b>Accessories</b>			
Injection needle			
Illustrated IFUs			
Patient's safety card			
Patient's wristband			
Stickers			