

ophtha futur®

gas

sf6 | c2f6 | c3f8



High purity

Adjustable gas
concentration

Environmentally
conscious system

36 months shelf life

ophthafutur[®] gas

technical data-sheet

Product Name	ophthafutur [®] sf6	ophthafutur [®] c2f6	ophthafutur [®] c3f8
Active compound	sulfur hexafluoride	hexafluoroethane	octafluoropropane
Composition			
Sum Formula	SF ₆	C ₂ F ₆	C ₃ F ₈
Purity	> 4.5	> 5.0	> 4.0
Mass fraction	> 99.995%	> 99.999%	> 99.99%
Non expansive gas concentration	20%	16%	12%
Potential tamponade time	6 days	15 days	30 days
Retention time	10-14 days	30-35 days	55-65 days
Properties			
Pure Gas Volume	15 ml	12 ml	9 ml
Total Syringe Volume - Mixing Device	60 ml	60 ml	60 ml
Shelf life	36 months	36 months	36 months
Presentation			
Container	Glass Syringe	Glass Syringe	Glass Syringe
Mixing Device	Plastic Syringe with Filter	Plastic Syringe with Filter	Plastic Syringe with Filter
Accessories			
Injection needle			
Illustrated IFUs			
Patient's safety card			
Patient's wristband			
Stickers			