

ophtha futur® octa | deca



High patient and
product safety

Multi-step proven & tested
ultra-purification process

Safe, biocompatible, sterile,
endotoxin-free products

Novel high performance
polymer syringes (HPPS)

36 months shelf life



ophthafutur[®] octa | deca

technical data-sheet

Product Name

Active compound

ophthafutur[®] octa
perfluorooctane

ophthafutur[®] deca
perfluorodecalin

Physical and Chemical Parameters

Sum Formula	C ₈ F ₁₈	C ₁₀ F ₁₈
Refractive Index n _D ²⁰	1.27	1.31
Density (20 °C)	1.77 g/ml	1.93 g/ml
Vapor pressure (at 37°C)	8.8 kPa	3.4 kPa
Compound Surface Tension	14 mN/m	19 mN/m
Interfacial Tension against Water	54 mN/m	53 mN/m

Chemical Control

Fluorination level	100 %	100 %
H-value	≤ 10 ppm	≤ 10 ppm
Additionally controlled Impurities (GC/MS)		
• concentration of C-H components	< 10 ppm	< 10 ppm
• 1H-PFO content	< 1 ppm	n.a.
• Perfluorofurane content	< 0.4%	n.a.

Microbial control

Product Sterility ensured by	Sterile Filtration	Sterile Filtration
Outer Surface Sterility ensured by	Steam Sterilization	Steam Sterilization
Bacterial Endotoxins ¹	≤ 0.2 EU/ml	≤ 0.2 EU/ml

¹ LAL test acc. to EP 2.6.14/USP <85>

Biocompatibility acc. to DIN EN ISO 10993

Cytotoxicity ²	Non cytotoxic	Non cytotoxic
Sensitization ³	Non sensitizing	Non sensitizing
Irritation ⁴	Non irritating	Non irritating

² Growth inhibition test of serum extract acc. to DIN EN ISO 10993-5:2009

³ Murine Local Lymph Node Assay (LLNA) acc. to DIN EN ISO 10993-10:2010

⁴ Acute eye-irritation on rabbits acc. to DIN EN ISO 10993-10:2010 Appendix B

Cell Toxicity (%)

H-value (ppm)	0	8	14	36	72	87	92
	≤ 10	175	350	700	1,400	2,100	2,800

Menz D-H, Feltgen N, Menz H, et al. How to ward off retinal toxicity of perfluorooctane and other perfluorocarbon liquids?

Invest Ophthalmol Vis. Sci. 2018; 59: 4841-4846.

Presentation

	Filling Volume	Container	Sterile Barrier	Shelf Life
ophthafutur [®] deca	5 ml / 7 ml	Vial	Pouch	36 months
ophthafutur [®] deca	5 ml / 7 ml	Syringe	Pouch	36 months
ophthafutur [®] octa	5 ml / 7 ml	Vial	Pouch	36 months
ophthafutur [®] octa	5 ml / 7 ml	Syringe	Pouch	36 months